

**Contenitore PS sovrapponibile****Hünersdorff**

Contenitore PS sovrapponibile. In polistirolo. Robusto, qualità industriale. Impilabile (in lunghezza e trasversalmente). Con certificato RAL. In PS ad alto impatto, resistente agli shock. Colori adatti per alimenti. Con la parte sotto liscia, adatto quindi al trasporto.



Code	Description	Packaging
LLG09033121	Scatola trasparente, colore rosso , largh. 100mm, lungh. 165/140mm, h 75mm	1 pz
LLG09033123	Scatola trasparente, colore blu, largh. 100mm, lungh. 165/140mm, h 75mm	1 pz
LLG09033131	Scatola trasparente, colore rosso, largh. 145mm, lungh. 232/200mm, h 125mm	1 pz
LLG09033133	Scatola trasparente, colore blu, largh. 145mm, lungh. 232/200mm, h 125mm	1 pz
LLG09033141	Scatola trasparente, colore rosso, largh. 200mm, lungh. 336/300mm, h 150mm	1 pz
LLG09033143	Scatola trasparente, colore blu, largh. 200mm, lungh. 336/300mm, h 150mm	1 pz



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

