

**Accessori per sonopuls HD 2200, in vetro****Bandelin**

Accessori per sonopuls HD 2200, in vetro. Celle per trattamento a ultrasuoni diretto. In vetro.



Code	Description	Packaging
LLG09650223	Cella con camicia raffreddamento, Ø25mm, 70ml, per HD 2200/3200	1 pz
LLG09650224	Cella a flusso continuo c/camicia raffredd., Ø25mm, 70ml, per HD 2200/3200	1 pz
LLG09650235	Cella con camicia raffreddamento, Ø13mm, 15ml, per HD 2070/2200/3100/3200	1 pz
LLG09650236	Cella flusso continuo c/camicia raffr. Ø13mm 15ml, x HD2070/2200/3100/3200	1 pz
LLG09650260	Cella Rosette per sonotrodi Ø max. 25 mm, 390 ml, per HD 2200/3200/3400	1 pz
LLG09650261	Cella Rosette per sonotrodi Ø max.13mm, 110 ml, per HD 2070/2200/3100/3200	1 pz
LLG09650264	Cella Rosette per sonotrodi Ø max. 25 mm, 660 ml, per HD 2200/3200/3400	1 pz
LLG09650267	Cella Rosette per sonotrodi Ø max.13 mm, 25 ml, per HD 2070/2200/3100/3200	1 pz
LLG09650268	Cella Rosette per sonotrodi Ø max.13 mm, 40 ml, per HD 2070/2200/3100/3200	1 pz



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

