

**Bicchieri, vetro Borosilicato 3.3, forma alta****Scherf Präzision**

Bicchieri, vetro Borosilicato 3.3, forma alta. DIN 12 331. ISO 3819. Con graduazioni approssimative e becco.



Code	Description	Packaging
LLG06270117	Bicchiera vetro borosilicato, forma alta, 2000 ml, Ø est. 120 mm, h 240 mm	10 pz
LLG06313628	Bicchiera, vetro borosilicato, forma alta, 3000 ml, Ø est. 135 mm, h 280 mm	1 pz
LLG09013920	Bicchiera, vetro borosilicato, forma alta, 100 ml, Ø est. 48 mm, h 80 mm	10 pz
LLG09013921	Bicchiera, vetro borosilicato, forma alta, 150 ml, Ø est. 54 mm, h 95 mm	10 pz
LLG09013922	Bicchiera, vetro borosilicato, forma alta, 250 ml, Ø est. 60 mm, h 120 mm	10 pz
LLG09013923	Bicchiera, vetro borosilicato, forma alta, 400 ml, Ø est. 70 mm, h 130 mm	10 pz
LLG09013924	Bicchiera, vetro borosilicato, forma alta, 25 ml, Ø est. 30 mm, h 55 mm	10 pz
LLG09013925	Bicchiera, vetro borosilicato, forma alta, 50 ml, Ø est. 38 mm, h 70 mm	10 pz
LLG09013930	Bicchiera, vetro borosilicato, forma alta, 600 ml, Ø est. 80 mm, h 150 mm	10 pz
LLG09013931	Bicchiera, vetro borosilicato, forma alta, 800 ml, Ø est. 90 mm, h 175 mm	10 pz
LLG09013932	Bicchiera, vetro borosilicato, forma alta, 1 L, Ø est. 95 mm, h 180 mm	10 pz

