


**Beute Erlenmeyer, vetro borosilicato 3.3, bocca stretta****Scherf Präzision**

Beute Erlenmeyer, vetro borosilicato 3.3, bocca stretta. Conforme DIN ISO 1773.
Con graduazioni e bordo svasato.



Code	Description	Packaging
LLG06270089	Beuta bocca stretta, 2000 ml, Ø collo 50 mm, Ø est. 166 mm, h 280 mm	10 pz
LLG06270118	Beuta bocca stretta, 3000 ml, Ø collo 50 mm, Ø est. 187 mm, h 310 mm	1 pz
LLG06270119	Beuta bocca stretta, 5000 ml, Ø collo 50 mm, Ø est. 220 mm, h 365 mm	1 pz
LLG09141079	Beuta bocca stretta, 25 ml, Ø collo 22 mm, Ø est. 42 mm, h 70 mm	10 pz
LLG09141080	Beuta bocca stretta, 50 ml, Ø collo 22 mm, Ø est. 51 mm, h 85 mm	10 pz
LLG09141081	Beuta bocca stretta, 100 ml Ø collo 22 mm, Ø est. 64 mm, h 105 mm	10 pz
LLG09141082	Beuta bocca stretta, 250 ml Ø collo 34 mm, Ø est. 85 mm, h 140 mm	10 pz
LLG09141083	Beuta bocca stretta, 500 ml Ø collo 34 mm, Ø est. 105 mm, h 175 mm	10 pz
LLG09141084	Beuta bocca stretta, 1000 ml, Ø collo 42 mm, Ø est. 131 mm, h 220 mm	10 pz
LLG09141085	Beuta bocca stretta, 200 ml, Ø collo 34 mm, Ø est. 79 mm, h 131 mm	10 pz
LLG09141086	Beuta bocca stretta, 300 ml, Ø collo 34 mm, Ø est. 87 mm, h 156 mm	10 pz



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

