


**Beute Erlenmeyer, vetro borosilicato 3.3, bocca stretta****Scherf Präzision**

Beute Erlenmeyer, vetro borosilicato 3.3, bocca stretta. Conforme DIN ISO 1773.  
Con graduazioni e bordo svasato.



| Code        | Description                                                          | Packaging |
|-------------|----------------------------------------------------------------------|-----------|
| LLG06270089 | Beuta bocca stretta, 2000 ml, Ø collo 50 mm, Ø est. 166 mm, h 280 mm | 10 pz     |
| LLG06270118 | Beuta bocca stretta, 3000 ml, Ø collo 50 mm, Ø est. 187 mm, h 310 mm | 1 pz      |
| LLG06270119 | Beuta bocca stretta, 5000 ml, Ø collo 50 mm, Ø est. 220 mm, h 365 mm | 1 pz      |
| LLG09141079 | Beuta bocca stretta, 25 ml, Ø collo 22 mm, Ø est. 42 mm, h 70 mm     | 10 pz     |
| LLG09141080 | Beuta bocca stretta, 50 ml, Ø collo 22 mm, Ø est. 51 mm, h 85 mm     | 10 pz     |
| LLG09141081 | Beuta bocca stretta, 100 ml Ø collo 22 mm, Ø est. 64 mm, h 105 mm    | 10 pz     |
| LLG09141082 | Beuta bocca stretta, 250 ml Ø collo 34 mm, Ø est. 85 mm, h 140 mm    | 10 pz     |
| LLG09141083 | Beuta bocca stretta, 500 ml Ø collo 34 mm, Ø est. 105 mm, h 175 mm   | 10 pz     |
| LLG09141084 | Beuta bocca stretta, 1000 ml, Ø collo 42 mm, Ø est. 131 mm, h 220 mm | 10 pz     |
| LLG09141085 | Beuta bocca stretta, 200 ml, Ø collo 34 mm, Ø est. 79 mm, h 131 mm   | 10 pz     |
| LLG09141086 | Beuta bocca stretta, 300 ml, Ø collo 34 mm, Ø est. 87 mm, h 156 mm   | 10 pz     |