


**Beute Erlenmeyer, vetro Borosilicato 3.3, bocca larga****Scherf Präzision**

Beute Erlenmeyer, vetro Borosilicato 3.3, bocca larga. Conforme DIN ISO 24450. Con graduazioni e bordo svasato.



Code	Description	Packaging
LLG06270181	Beuta bocca larga, 2000 ml, Ø collo 72 mm, Ø est. 153 mm, h 280 mm	10 pz
LLG09141180	Beuta bocca larga, 50 ml, Ø collo 34 mm, Ø est. 51 mm, h 85 mm	10 pz
LLG09141181	Beuta bocca larga, 100 ml, Ø collo 34 mm, Ø est. 64 mm, h 105 mm	10 pz
LLG09141182	Beuta bocca larga, 250 ml, Ø collo 50 mm, Ø est. 85 mm, h 135 mm	10 pz
LLG09141183	Beuta bocca larga, 300 ml, Ø collo 50 mm, Ø est. 87 mm, h 156 mm	10 pz
LLG09141184	Beuta bocca larga, 500 ml, Ø collo 50 mm, Ø est. 105 mm, h 170 mm	10 pz
LLG09141185	Beuta bocca larga, 1000 ml, Ø collo 50 mm, Ø est. 131 mm, h 215 mm	10 pz
LLG09141186	Beuta bocca larga, 25 ml, Ø collo 32 mm, Ø est. 42 mm, h 70 mm	10 pz
LLG09141187	Beuta bocca larga, 200 ml, Ø collo 50 mm, Ø est. 79 mm, h 131 mm	10 pz