

**Tubo in PTFE****Deutsch & Neumann**

Tubo, PTFE.

Politetrafluoretilene (PTFE), bianco, senza pressione, resistente al calore fino a 270 °C.



Code	Description	Packaging
LLG06070190	Tubo, in PTFE, Ø int. 1,00 mm, Ø est. 3,0 mm, spessore pareti 1,00 mm ATTENZIONE! ORDINE MINIMO 5 CONF. O MULTIPLI DI 5 CONF.	1 pz



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

