

**Bacinelle per laboratorio in PS****Kartell**

Vassoi e vaschette, PS.

Adatte per usi alimentari ed usi generali di laboratorio. Polistirene antiurto.



Code	Description	Packaging
LLG06205063	Vassoi e vaschette, PS antiurto, dimensioni esterne (LxPxH) 408 x 300 x 21 mm, dimensioni interne (LxPxH) 380 x 270 x 10 mm	1 pz
LLG06206152	Vassoi e vaschette, PS antiurto, dimensioni esterne (LxPxH) 254 x 254 x 21 mm, dimensioni interne (LxPxH) 220 x 220 x 10 mm	1 pz
LLG06224951	Vassoi e vaschette, PS antiurto, dimensioni esterne (LxPxH) 409 x 300 x 42 mm, dimensioni interne (LxPxH) 360 x 250 x 30 mm	1 pz
LLG06224952	Vassoi e vaschette, PS antiurto, dimensioni esterne (LxPxH) 353 x 252 x 21 mm, dimensioni interne (LxPxH) 320 x 220 x 10 mm	1 pz.
LLG06224986	Vassoi e vaschette, PS antiurto, dimensioni esterne (LxPxH) 201 x 151 x 41 mm, dimensioni interne (LxPxH) 150 x 100 x 30 mm	1 pz
LLG06225054	Vassoi e vaschette, PS antiurto, dimensioni esterne (LxPxH) 303 x 151 x 42 mm, dimensioni interne (LxPxH) 260 x 110 x 30 mm	1 pz
LLG06225055	Vassoi e vaschette, PS antiurto, dimensioni esterne (LxPxH) 353 x 253 x 41 mm, dimensioni interne (LxPxH) 310 x 210 x 30 mm	1 pz
LLG06225141	Vassoi e vaschette, PS antiurto, dimensioni esterne (LxPxH) 303 x 151 x 21 mm, dimensioni interne (LxPxH) 270 x 120 x 10 mm	1 pz
LLG06227810	Vassoi e vaschette, PS antiurto, dimensioni esterne (LxPxH) 352 x 252 x 81 mm, dimensioni interne (LxPxH) 290 x 190 x 70 mm	1 pz
LLG06227811	Vassoi e vaschette, PS antiurto, dimensioni esterne (LxPxH) 408 x 299 x 81 mm, dimensioni interne (LxPxH) 340 x 240 x 70 mm	1 pz
LLG06228312	Vassoi e vaschette, PS antiurto, dimensioni esterne (LxPxH) 202 x 151 x 81 mm, dimensioni interne (LxPxH) 140 x 90 x 70 mm	1 pz
LLG06233344	Vassoi e vaschette, PS antiurto, dimensioni esterne (LxPxH) 201 x 151 x 21 mm, dimensioni interne (LxPxH) 170 x 110 x 10 mm	1 pz





CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

