

**Bottiglia forma quadrata bocca stretta serie 310, PE-HD, PVC****Kautex**

Bottiglia forma quadrata bocca stretta serie 310, PE-HD, PVC. Bottiglia quadrata, in PE-HD, naturale, translucida. Bottiglia rigida quadrata, in PVC, trasparente. Si prega ordinare a parte il tappo a vite.



Code	Description	Packaging
LLG09072808	Bottiglia in PE-HD, 100 ml, 40,5x32,0x110,5 h mm, Ø filettatura est. 18 mm	1 pz
LLG09072814	Bottiglia in PE-HD, 250 ml, 53,5x42,0x146,5 h mm, Ø filettatura est. 25 mm	1 pz
LLG09072819	Bottiglia in PE-HD, 500 ml, 63,5x54,5x180 h mm, Ø filettatura est. 25 mm	1 pz
LLG09072824	Bottiglia in PE-HD, 1000 ml, 84,0x66x222,5 h mm, Ø filettatura est. 25 mm	1 pz
LLG09072914	Bottiglia in PVC, 250 ml, 59x38x151,5 h mm, Ø filettatura est. 25 mm	1 pz
LLG09072919	Bottiglia in PVC, 500 ml, 70x50x189 h mm, Ø filettatura est. 25 mm	1 pz
LLG09072924	Bottiglia in PVC, 1000 ml, 87,5x61,0x233 h mm, Ø filettatura est. 25 mm	1 pz



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

