

**Bottiglia forma quadrata bocca stretta serie 310, PE-HD, PVC****Kautex**

Bottiglia forma quadrata bocca stretta serie 310, PE-HD, PVC. Bottiglia quadrata, in PE-HD, naturale, translucida. Bottiglia rigida quadrata, in PVC, trasparente. Si prega ordinare a parte il tappo a vite.



| Code        | Description                                                                | Packaging |
|-------------|----------------------------------------------------------------------------|-----------|
| LLG09072808 | Bottiglia in PE-HD, 100 ml, 40,5x32,0x110,5 h mm, Ø filettatura est. 18 mm | 1 pz      |
| LLG09072814 | Bottiglia in PE-HD, 250 ml, 53,5x42,0x146,5 h mm, Ø filettatura est. 25 mm | 1 pz      |
| LLG09072819 | Bottiglia in PE-HD, 500 ml, 63,5x54,5x180 h mm, Ø filettatura est. 25 mm   | 1 pz      |
| LLG09072824 | Bottiglia in PE-HD, 1000 ml, 84,0x66x222,5 h mm, Ø filettatura est. 25 mm  | 1 pz      |
| LLG09072914 | Bottiglia in PVC, 250 ml, 59x38x151,5 h mm, Ø filettatura est. 25 mm       | 1 pz      |
| LLG09072919 | Bottiglia in PVC, 500 ml, 70x50x189 h mm, Ø filettatura est. 25 mm         | 1 pz      |
| LLG09072924 | Bottiglia in PVC, 1000 ml, 87,5x61,0x233 h mm, Ø filettatura est. 25 mm    | 1 pz      |



**CARLO ERBA Reagents S.r.l.**  
 VIA R. MERENDI, 22  
 20007 CORNAREDO (MI)  
 TEL. +39 02 93 99 190  
 FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

