

**Elettrodi pH, digital IDS SenTix 950****WTW**

Elettrodi pH, digital IDS SenTix® 950. Caratteristiche: Range temperatura: 0 ... 80 °C - Range pH: 0 ... 14 pH - Lunghezza / diam. : 120mm / 12mm. Resistenza membrana a 25°C: 1GW.



Code	Description	Packaging
LLG07656698	Elettrodi pH, digital IDS SenTix 950	1 pz



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

