



Test in cuvetta per analisi dell'acqua

HACH LANGE

Test in cuvetta per analisi dell'acqua



Code	Description	Packaging
LLG06231183	Test in cuvetta ammonio 2,0-47,0 mg/l NH ₄ -N, 25 test	1 pz
LLG06231184	Cuvette Tipo LCK349, per fosfato totale, 0.05 - 1.50mg/L PO ₄ -P	1 pz
LLG06231242	Test in cuvetta per azoto totale 20-100 mg/L TNb, 25 test	1 pz
LLG06231248	Test in cuvetta per cloruro, 1 - 70 mg/l / 70 - 1000 mg/l Cl, 24 test	1 pz
LLG06242605	Test in cuvetta per fosfato (orto/totale), 2,0 - 20,0 mg/l PO ₄ -P, 25 test	1 pz
LLG06802583	Test in cuvetta per COD, 100 - 2000 mg/l O ₂ , 25 test	1 pz
LLG06803024	Tensioattivi anionici, test in cuvetta, 0,05-2,0 mg/L	1 pz



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

