

**Tappo in silicone****Deutsch & Neumann**

Tappo in silicone. Prodotti con tubo DURAN®, con rubinetto a valvola per dosaggio, con oliva in vetro (angolo 90°), con cono NS 14/23, lunghezza tubo dalla cima al giunto: circa 150mm.



Code	Description	Packaging
LLG09230503	Tappo in silicone Ø base 3,5 mm, Ø sopra 6,5 mm, h 15 mm	1 pz
LLG09230505	Tappo in silicone Ø base 5,0 mm, Ø sopra 9,0 mm, h 20 mm	1 pz
LLG09230508	Tappo in silicone Ø base 8,0 mm, Ø sopra 12,0 mm, h 20 mm	1 pz
LLG09230510	Tappo in silicone Ø base 10,5 mm, Ø sopra 14,5 mm, h 20 mm	1 pz
LLG09230512	Tappo in silicone Ø base 12,5 mm, Ø sopra 16,5 mm, h 20 mm	1 pz
LLG09230514	Tappo in silicone Ø base 14,0 mm, Ø sopra 18,0 mm, h 20 mm	1 pz
LLG09230517	Tappo in silicone Ø base 17,0 mm, Ø sopra 22,0 mm, h 25 mm	1 pz
LLG09230518	Tappo in silicone Ø base 18,0 mm, Ø sopra 24,0 mm, h 30 mm	1 pz
LLG09230521	Tappo in silicone Ø base 21,0 mm, Ø sopra 27,0 mm, h 30 mm	1 pz
LLG09230523	Tappo in silicone Ø base 23,0 mm, Ø sopra 29,0 mm, h 30 mm	1 pz
LLG09230526	Tappo in silicone Ø base 26,0 mm, Ø sopra 32,0 mm, h 30 mm	1 pz
LLG09230529	Tappo in silicone Ø base 29,0 mm, Ø sopra 35,0 mm, h 30 mm	1 pz
LLG09230531	Tappo in silicone Ø base 31,0 mm, Ø sopra 38,0 mm, h 35 mm	1 pz
LLG09230536	Tappo in silicone Ø base 36,0 mm, Ø sopra 44,0 mm, h 40 mm	1 pz
LLG09230541	Tappo in silicone Ø base 41,0 mm, Ø sopra 49,0 mm, h 40 mm	1 pz
LLG09230547	Tappo in silicone Ø base 47,0 mm, Ø sopra 55,0 mm, h 40 mm	1 pz
LLG09230550	Tappo in silicone Ø base 50,5 mm, Ø sopra 59,5 mm, h 45 mm	1 pz
LLG09230556	Tappo in silicone Ø base 56,0 mm, Ø sopra 65,0 mm, h 45 mm	1 pz
LLG09230560	Tappo in silicone Ø base 60,0 mm, Ø sopra 70,0 mm, h 50 mm	1 pz
LLG09230564	Tappo in silicone Ø base 64,5 mm, Ø sopra 75,5 mm, h 55 mm	1 pz



Code	Description	Packaging
LLG09230571	Tappo in silicone Ø base 71,0 mm, Ø sopra 83,0 mm, h 60 mm	1 pz
LLG09230579	Tappo in silicone Ø base 79,0 mm Ø sopra 92,0 mm h 65 mm	1 pz



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

