

**Tappo in silicone****Deutsch & Neumann**

Tappo in silicone. Prodotti con tubo DURAN®, con rubinetto a valvola per dosaggio, con oliva in vetro (angolo 90°), con cono NS 14/23, lunghezza tubo dalla cima al giunto: circa 150mm.



| Code        | Description                                                | Packaging |
|-------------|------------------------------------------------------------|-----------|
| LLG09230503 | Tappo in silicone Ø base 3,5 mm, Ø sopra 6,5 mm, h 15 mm   | 1 pz      |
| LLG09230505 | Tappo in silicone Ø base 5,0 mm, Ø sopra 9,0 mm, h 20 mm   | 1 pz      |
| LLG09230508 | Tappo in silicone Ø base 8,0 mm, Ø sopra 12,0 mm, h 20 mm  | 1 pz      |
| LLG09230510 | Tappo in silicone Ø base 10,5 mm, Ø sopra 14,5 mm, h 20 mm | 1 pz      |
| LLG09230512 | Tappo in silicone Ø base 12,5 mm, Ø sopra 16,5 mm, h 20 mm | 1 pz      |
| LLG09230514 | Tappo in silicone Ø base 14,0 mm, Ø sopra 18,0 mm, h 20 mm | 1 pz      |
| LLG09230517 | Tappo in silicone Ø base 17,0 mm, Ø sopra 22,0 mm, h 25 mm | 1 pz      |
| LLG09230518 | Tappo in silicone Ø base 18,0 mm, Ø sopra 24,0 mm, h 30 mm | 1 pz      |
| LLG09230521 | Tappo in silicone Ø base 21,0 mm, Ø sopra 27,0 mm, h 30 mm | 1 pz      |
| LLG09230523 | Tappo in silicone Ø base 23,0 mm, Ø sopra 29,0 mm, h 30 mm | 1 pz      |
| LLG09230526 | Tappo in silicone Ø base 26,0 mm, Ø sopra 32,0 mm, h 30 mm | 1 pz      |
| LLG09230529 | Tappo in silicone Ø base 29,0 mm, Ø sopra 35,0 mm, h 30 mm | 1 pz      |
| LLG09230531 | Tappo in silicone Ø base 31,0 mm, Ø sopra 38,0 mm, h 35 mm | 1 pz      |
| LLG09230536 | Tappo in silicone Ø base 36,0 mm, Ø sopra 44,0 mm, h 40 mm | 1 pz      |
| LLG09230541 | Tappo in silicone Ø base 41,0 mm, Ø sopra 49,0 mm, h 40 mm | 1 pz      |
| LLG09230547 | Tappo in silicone Ø base 47,0 mm, Ø sopra 55,0 mm, h 40 mm | 1 pz      |
| LLG09230550 | Tappo in silicone Ø base 50,5 mm, Ø sopra 59,5 mm, h 45 mm | 1 pz      |
| LLG09230556 | Tappo in silicone Ø base 56,0 mm, Ø sopra 65,0 mm, h 45 mm | 1 pz      |
| LLG09230560 | Tappo in silicone Ø base 60,0 mm, Ø sopra 70,0 mm, h 50 mm | 1 pz      |
| LLG09230564 | Tappo in silicone Ø base 64,5 mm, Ø sopra 75,5 mm, h 55 mm | 1 pz      |



| Code        | Description                                                | Packaging |
|-------------|------------------------------------------------------------|-----------|
| LLG09230571 | Tappo in silicone Ø base 71,0 mm, Ø sopra 83,0 mm, h 60 mm | 1 pz      |
| LLG09230579 | Tappo in silicone Ø base 79,0 mm Ø sopra 92,0 mm h 65 mm   | 1 pz      |



**CARLO ERBA Reagents S.r.l.**  
 VIA R. MERENDI, 22  
 20007 CORNAREDO (MI)  
 TEL. +39 02 93 99 190  
 FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

