



Accessori per turbidimetro AL250T-IR / AL450T-IR

Tintometer

Accessori per turbidimetro AL250T-IR / AL450T-IR.



Code	Description	Packaging
LLG06223191	Set di 12 cuvette portacampione vuote, Ø 24 mm	1 pz
LLG09920493	Set di standard di turbidità T-CAL (<0,1, 20, 200, 800 NTU)	1 pz



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

