

**Bicchieri di misura, PP****LLG Labware**

Bicchieri di misura, PP. Conformi ISO 7056. Graduazioni stampate blu.



| Code        | Description                                                        | Packaging |
|-------------|--------------------------------------------------------------------|-----------|
| LLG09013670 | Bicchieri di Griffin, PP, 25 ml, div. 5 ml, Ø 32 mm, h 47 mm       | 2 pz      |
| LLG09013671 | Bicchieri di Griffin, PP, 50 ml, div. 10 ml, Ø 39 mm, h 60 mm      | 2 pz      |
| LLG09013672 | Bicchieri di Griffin, PP, 100 ml, div. 10 ml, Ø 50 mm, h 70 mm     | 2 pz      |
| LLG09013673 | Bicchieri di Griffin, PP, 250 ml, div. 25 ml, Ø 67 mm, h 93 mm     | 2 pz      |
| LLG09013674 | Bicchieri di Griffin, PP, 500 ml, div. 50 ml, Ø 88 mm, h 119 mm    | 2 pz      |
| LLG09013675 | Bicchieri di Griffin, PP, 1000 ml, div. 50 ml, Ø 109 mm, h 147 mm  | 2 pz      |
| LLG09013676 | Bicchieri di Griffin, PP, 2000 ml, div. 100 ml, Ø 133 mm, h 184 mm | 1 pz      |



**CARLO ERBA Reagents S.r.l.**  
 VIA R. MERENDI, 22  
 20007 CORNAREDO (MI)  
 TEL. +39 02 93 99 190  
 FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

