

**Bicchieri di misura, PP****LLG Labware**

Bicchieri di misura, PP. Conformi ISO 7056. Graduazioni stampate blu.



Code	Description	Packaging
LLG09013670	Bicchiere di Griffin, PP, 25 ml, div. 5 ml, Ø 32 mm, h 47 mm	2 pz
LLG09013671	Bicchiere di Griffin, PP, 50 ml, div. 10 ml, Ø 39 mm, h 60 mm	2 pz
LLG09013672	Bicchiere di Griffin, PP, 100 ml, div. 10 ml, Ø 50 mm, h 70 mm	2 pz
LLG09013673	Bicchiere di Griffin, PP, 250 ml, div. 25 ml, Ø 67 mm, h 93 mm	2 pz
LLG09013674	Bicchiere di Griffin, PP, 500 ml, div. 50 ml, Ø 88 mm, h 119 mm	2 pz
LLG09013675	Bicchiere di Griffin, PP, 1000 ml, div. 50 ml, Ø 109 mm, h 147 mm	2 pz
LLG09013676	Bicchiere di Griffin, PP, 2000 ml, div. 100 ml, Ø 133 mm, h 184 mm	1 pz

