

**Provette, vetro soda-lime AR®****LLG Labware**

Provette con fondo tondo in vetro AR®. Resistenza idrolitica: Classe HGB 3 secondo ISO 719, Tipo III conforme a C- USP.



Code	Description	Packaging
LLG06201034	Provette vetro soda-lime AR, Ø 8 mm, h 70 mm, con bordo, spess. 0,5 mm	100 pz
LLG06267806	Provetta in vetro, Ø30x100 h mm, spess. 1 mm, senza bordo	100 pz
LLG06800382	Provette vetro soda-lime AR, Ø 16 mm, h 130 mm, con bordo, spess. 0,6 mm	100 pz
LLG06800620	Provette vetro soda-lime AR, Ø 14 mm, h 130 mm, con bordo, spess. 0,5 mm	100 pz
LLG06800621	Provette vetro soda-lime AR, Ø 25 mm, h 150 mm, con bordo, spess. 0,9 mm	50 pz
LLG06801073	Provette vetro soda-lime AR, Ø 10 mm, h 100 mm, con bordo, spess. 0,5 mm	100 pz
LLG07600324	Provette vetro soda-lime AR, Ø 30 mm, h 200 mm, con bordo, spess. 0,9 mm	50 pz
LLG07600405	Provette vetro soda-lime AR, Ø 16 mm, h 160 mm, con bordo, spess. 0,6 mm	100 pz
LLG07607198	Provette vetro soda-lime AR, Ø 12 mm, h 100 mm, con bordo, spess. 0,5 mm	100 pz
LLG07607778	Provette vetro soda-lime AR, Ø 20 mm, h 180 mm, con bordo, spess. 0,6 mm	100 pz
LLG07651871	Provette vetro soda-lime AR, Ø 16 mm, h 120 mm, senza bordo, spess. 1,0 mm	100 pz
LLG07651872	Provette vetro soda-lime AR, Ø 15 mm, h 100 mm, senza bordo, spess. 1,0 mm	100 pz
LLG07652552	Provette vetro soda-lime AR, Ø 18 mm, h 130 mm, senza bordo spess. 1,0 mm	100 pz



