




Accessori per Bagni a secco AccuBlock

Corning

Accessori per Bagni a secco AccuBlock#.



Code	Description	Packaging
LLG06224191	Blocco per 24 provette da 1,5 ml	1 pz
LLG06233664	Blocco per 20 provette da 2,0 ml	1 pz
LLG06235417	Blocco per 24 provette da 0,5 ml	1 pz
LLG06239530	Alublock per provette 12x15/16 mm	1 pz
LLG06241353	Alublock per provette 20x12 mm	1 pz
LLG06242290	Alublock per provette 20x10 mm	1 pz
LLG06251238	Alublock per provette 20x13 mm	1 pz
LLG07627809	Blocco riscaldante	1 pz
LLG09945763	Blocco per 48 provette PCR 0,2ml o strip 6 x 0,2ml	1 pz
LLG09945764	Blocco per 12 x 15 ml provette per centrifuga	1 pz
LLG09945765	Blocco per 5 x 50 ml provette per centrifuga	1 pz

