



Bottiglia a bocca stretta, PE-LD

LLG Labware

Bottiglia a bocca stretta, PE-LD. Trasparente, con tappo.



Code	Description	Packaging
LLG09073500	Bottiglia collo stretto, 50 ml Ø 37 mm, h 78 mm, Ø filettatura est. 18 mm	50 pz
LLG09073501	Bottiglia collo stretto, 100 ml Ø 47 mm, h 105 mm, Ø filettatura est. 18 mm	20 pz
LLG09073502	Bottiglia collo stretto, 250 ml Ø 63 mm, h 119 mm, Ø filettatura est. 25 mm	10 pz
LLG09073503	Bottiglia collo stretto, 500 ml Ø 78 mm, h 152 mm, Ø filettatura est. 25 mm	10 pz
LLG09073504	Bottiglia collo stretto, 1 L Ø 96 mm, h 206 mm, Ø filettatura est. 28 mm	4 pz
LLG09073505	Bottiglia collo stretto, 2 L Ø 120 mm, h 264 mm, Ø filettatura est. 28 mm	1 pz



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

