



Vaschetta tipo HELLENDahl

Glaswarenfabrik Karl Hecht

Vaschetta tipo HELLendhal. Per 16 vetrini portaoggetto da 26 x 76 mm. In vetro, con coperchio.



Code	Description	Packaging
LLG09163277	VASCHETTA tipo HELLEND AHL, Largh. 35 mm, Lungh. 85 mm, Altezza 95 mm	1 pz
LLG09163278	VASCHETTA tipo HELLEND AHL, Largh. 55 mm, Lungh. 60 mm, Altezza 105 mm	1 pz



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

