



Pinzette per uso generale, acciaio 18/10

LLG Labware

Pinzette per uso generale, acciaio 18/10. Con pinze a molla.



Code	Description	Packaging
LLG04008471	Pinzette in acciaio, lungh. 105 mm, smussate - dritte	1 pz
LLG04008472	Pinzette in acciaio, lungh. 115 mm, smussate - dritte	1 pz
LLG04008473	Pinzette in acciaio, lungh. 130 mm, smussate - dritte	1 pz
LLG04008474	Pinzette in acciaio, lungh. 160 mm, smussate - dritte	1 pz
LLG04008475	Pinzette in acciaio, lungh. 200 mm, smussate - dritte	1 pz
LLG04008476	Pinzette in acciaio, lungh. 250 mm, smussate - dritte	1 pz
LLG04008477	Pinzette in acciaio, lungh. 300 mm, smussate - dritte	1 pz
LLG04008478	Pinzette in acciaio, lungh. 115 mm, smussate - dritte, dentellate	1 pz
LLG04008479	Pinzette in acciaio, lungh. 130 mm, smussate - dritte, dentellate	1 pz
LLG04008480	Pinzette in acciaio, lungh. 145 mm, smussate - dritte, dentellate	1 pz
LLG04008481	Pinzette in acciaio, lungh. 160 mm, smussate - dritte, dentellate	1 pz
LLG06255826	Pinzette in acciaio, lungh. 145 mm, smussate - dritte	1 pz