



## LLG-Combustion boats, porcelain, with eyelet

LLG Labware

### LLG-Combustion boats, porcelain, with eyelet

Unglazed. With eyelet for an easy removal out the furnace after the annealing process.



| Code        | Description                             | Packaging |
|-------------|-----------------------------------------|-----------|
| LLG06243229 | LLG-Melting boats with eyelet 214/890x1 | 1 pz      |
| LLG09250940 | LLG-Melting boats,unglazed,with eyelet  | 1 pz      |
| LLG09250941 | LLG-Melting boats,unglazed,with eyelet  | 1 pz      |
| LLG09250942 | LLG-Melting boats,unglazed,with eyelet  | 1 pz      |
| LLG09250943 | LLG-Melting boats,unglazed,with eyelet  | 1 pz      |



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

