



**Filtri in cellulosa per ASE 100/300**

**Restek**

Filtri per ASE 100/300 (solo 34 mL, 66 mL e 100 mL), in cellulosa, 30 mm.



| Code        | Description                         | Packaging |
|-------------|-------------------------------------|-----------|
| LLG06496550 | Filtri in cellulosa per ASE 100/300 | 100 pz    |



**CARLO ERBA Reagents S.r.l.**  
 VIA R. MERENDI, 22  
 20007 CORNAREDO (MI)  
 TEL. +39 02 93 99 190  
 FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

