

**Tubi Nessler, senza becco, vetro borosilicato 3.3****NORMAX**

Tubi Nessler, senza becco, vetro borosilicato 3.3. Ad esempio per uso con il sistema Nessler di Lovibond®.



Code	Description	Packaging
LLG09170621	Tubo Nessler, forma alta, graduato a 50 ml, cap. 50 ml, (Øxh) 20x300 mm	1 pz
LLG09170622	Tubo Nessler, forma alta, graduato a 100 ml, cap. 100 ml, (Øxh) 24x375 mm	1 pz
LLG09170623	Tubo Nessler, forma bassa, graduato a 50 ml, cap. 50 ml, (Øxh) 25x175 mm	1 pz
LLG09170624	Tubo Nessler, forma bassa, graduato a 100 ml, cap. 100 ml, (Øxh) 32x200 mm	1 pz



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

