

**Tubi Nessler, senza becco, vetro borosilicato 3.3****NORMAX**

Tubi Nessler, senza becco, vetro borosilicato 3.3. Ad esempio per uso con il sistema Nessler di Lovibond®.



Code	Description	Packaging
LLG09170621	Tubo Nessler, forma alta, graduato a 50 ml, cap. 50 ml, (Øxh) 20x300 mm	1 pz
LLG09170622	Tubo Nessler, forma alta, graduato a 100 ml, cap. 100 ml, (Øxh) 24x375 mm	1 pz
LLG09170623	Tubo Nessler, forma bassa, graduato a 50 ml, cap. 50 ml, (Øxh) 25x175 mm	1 pz
LLG09170624	Tubo Nessler, forma bassa, graduato a 100 ml, cap. 100 ml, (Øxh) 32x200 mm	1 pz



**CARLO ERBA Reagents S.r.l.**  
 VIA R. MERENDI, 22  
 20007 CORNAREDO (MI)  
 TEL. +39 02 93 99 190  
 FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

