



**Sessole, HDPE**

**LLG Labware**

Sessole, HDPE. Per uso industriale. Adatte per alimenti.



| Code        | Description                                    | Packaging |
|-------------|------------------------------------------------|-----------|
| LLG04682195 | Sessole, HDPE, 1000 ml, Lungh. 350 mm          | 5 pz      |
| LLG04682196 | Sessole, HDPE, 1500 ml, Lungh. 400 mm          | 5 pz      |
| LLG04682197 | Sessole, HDPE, Capacità 250 ml, Lungh. 250 mm  | 5 pz      |
| LLG04682198 | Sessole, HDPE, 500 ml, Lungh. 300 mm           | 5 pz      |
| LLG06240471 | Sessole, HDPE, lungh. 350 mm, cap. 1000 ml     | 1 pz      |
| LLG06321690 | Sessole, HDPE, Capacità 1000 ml, Lungh. 350 mm | 25pz      |
| LLG06321691 | Sessole, HDPE, Capacità 250 ml, Lungh. 250 mm  | 25pz      |
| LLG06321692 | Sessole, HDPE, Capacità 500 ml, Lungh. 300 mm  | 25pz      |
| LLG09201091 | Sessole, HDPE, lungh. 250 mm, cap. 250 ml      | 1 pz      |
| LLG09201092 | Sessole, HDPE, lungh. 300 mm, cap. 500 ml      | 1 pz      |
| LLG09201093 | Sessole, HDPE, lungh. 400 mm, cap. 1500 ml     | 1 pz      |