



Sessole, HDPE

LLG Labware

Sessole, HDPE. Per uso industriale. Adatte per alimenti.



Code	Description	Packaging
LLG04682195	Sessole, HDPE, 1000 ml, Lungh. 350 mm	5 pz
LLG04682196	Sessole, HDPE, 1500 ml, Lungh. 400 mm	5 pz
LLG04682197	Sessole, HDPE, Capacità 250 ml, Lungh. 250 mm	5 pz
LLG04682198	Sessole, HDPE, 500 ml, Lungh. 300 mm	5 pz
LLG06240471	Sessole, HDPE, lungh. 350 mm, cap. 1000 ml	1 pz
LLG06321690	Sessole, HDPE, Capacità 1000 ml, Lungh. 350 mm	25pz
LLG06321691	Sessole, HDPE, Capacità 250 ml, Lungh. 250 mm	25pz
LLG06321692	Sessole, HDPE, Capacità 500 ml, Lungh. 300 mm	25pz
LLG09201091	Sessole, HDPE, lungh. 250 mm, cap. 250 ml	1 pz
LLG09201092	Sessole, HDPE, lungh. 300 mm, cap. 500 ml	1 pz
LLG09201093	Sessole, HDPE, lungh. 400 mm, cap. 1500 ml	1 pz

