



Colonne preparative VarioPrep NUCLEODUR C18 Gravity

Macherey Nagel

Colonne preparative VarioPrep NUCLEODUR® C18 Gravity
NUCLEODUR® C18 Gravity - C8 Gravity fase ad alta densità non polare

- disponibile come Octadecyl- (C18 - USP L1) e Octyl- (C8 - USP L7)
 - Pori 110 Å; dimensioni particelle 1,8 µm, 3 µm e 5 µm per C18, 1,8 e 5 µm per C8
 - particelle di 7, 10, 12 e 16 µm per separazioni preparative disponibili su richiesta.
 - Contenuto di carbonio 18 % C per C18, 11 % C per C8
 - ideale per sviluppo di metodi
 - permette HPLC a pH estremi (pH 1 - 11)
 - adatte per **LC/MS** grazie alle basse caratteristiche di legante
 - raccomandate per un vasto campo di applicazioni sofisticate.
 - Classi di prodotti separati velocemente: prodotti farmaceutici, es. analgesici, antiinfiammatori, antidepressivi; erbicidi; fitofarmaci; immunosoppressori.
- Dimensioni particelle 5 µm, 18% C



Code	Description	Packaging
MNL76209780	VP 10/8 NUCLEODUR C8 Gravity VarioPrep HPLC guard column cartridge VP 10/8 NUCLEODUR C8 Gravity length: 10 mm, ID: 8 mm	2 pz
MNL762103100	Colonne preparative VarioPrep NUCLEODUR C18 Gravity, Ø int. 10 mm, Lungh. 50 mm	1 pz
MNL762103210	VP 50/21 NUCLEODUR C18 Gravity, 5 µm VarioPrep HPLC column VP 50/21 NUCLEODUR C18 Gravity, 5 µm *** with adjustable Endfitting *** length: 50 mm, ID: 21 mm	1 pz
MNL762109100	VP 125/10 NUCLEODUR C18 Gravity, 5 µm VarioPrep HPLC column VP 125/10 NUCLEODUR C18 Gravity, 5 µm Length: 125 mm, ID: 10 mm	1 pz
MNL762109210	VP 125/21 NUCLEODUR C18 Gravity, 5 µm VarioPrep HPLC column VP 125/21 NUCLEODUR C18 Gravity, 5 µm Length: 125 mm, ID: 21 mm	1 pz
MNL762113100	Colonne preparative VarioPrep NUCLEODUR C18 Gravity, Ø int. 10 mm, Lungh. 250 mm	1 pz
MNL762113210	VP 250/21 NUCLEODUR C18 Gravity, 5 µm VarioPrep HPLC column VP 250/21 NUCLEODUR C18 Gravity, 5 µm *** with adjustable Endfitting *** length: 250 mm, ID: 21 mm	1 pz
MNL762160160	VP 10/16 NUCLEODUR C18 Gravity VarioPrep HPLC guard column cartridge (short) VP 10/16 NUCLEODUR C18 Gravity length: 10 mm, ID: 16 mm	2 pz
MNL76216080	Colonne preparative VarioPrep NUCLEODUR C18 Gravity, Ø int. 8* mm, Lungh. 10 mm	2 pz

