



Colonne capillari per analisi solventi e alcoli PERMABOND CW 20 M

Macherey Nagel

Colonne capillari per analisi solventi e alcoli PERMABOND CW 20 M.



Code	Description	Packaging
MNL72306010	Capillary column PERMABOND CW 20 M ID: 0.25 mm, film thickness: 0.25 µm length: 10 m	1 pz
MNL72306025	Capillary column PERMABOND CW 20 M ID: 0.25 mm, film thickness: 0.25 µm length: 25 m	1 pz
MNL72306030	Capillary column PERMABOND CW 20 M ID: 0.25 mm, film thickness: 0.25 µm length: 30 m	1 pz
MNL72306050	Capillary column PERMABOND CW 20 M ID: 0.25 mm, film thickness: 0.25 µm length: 50 m	1 pz
MNL72306060	Capillary Column PERMABOND CW 20 M ID: 0.25 mm, film thickness: 0.25 µm Length: 60 m	1 pz
MNL72306410	Capillary column PERMABOND CW 20 M ID: 0.10 mm, film thickness: 0.10 µm Length: 10 m	1 pz
MNL72329610	Capillary column PERMABOND CW 20 M ID: 0.32 mm, film thickness: 0.5 µm length: 10 m	1 pz
MNL72329625	Capillary column PERMABOND CW 20 M ID: 0.32 mm, film thickness: 0.5 µm length: 25 m	1 pz
MNL72329630	Capillary column PERMABOND CW 20 M ID: 0.32 mm, film thickness: 0.5 µm length: 30 m	1 pz
MNL72329650	Capillary column PERMABOND CW 20 M ID: 0.32 mm, film thickness: 0.5 µm length: 50 m	1 pz
MNL72329660	Capillary column PERMABOND CW 20 ID: 0.32 mm, film thickness: 0.5 µm length: 60 m	1 pz
MNL72332110	Capillary column PERMABOND CW 20 M ID: 0.32 mm, film thickness: 0.25 µm length: 10 m	1 pz
MNL72332125	Capillary column PERMABOND CW 20 M ID: 0.32 mm, film thickness: 0.25 µm length: 25 m	1 pz
MNL72332130	Capillary column PERMABOND CW 20 M ID: 0.32 mm, film thickness: 0.25 µm Length: 30 m	1 pz
MNL72332150	Capillary column PERMABOND CW 20 M ID: 0.32 mm, film thickness: 0.25 µm length: 50 m	1 pz
MNL72332160	Capillary column PERMABOND CW 20 M ID: 0.32 mm, film thickness: 0.25 µm length: 60 m Special manufacturing	1 pz
MNL72351510	Capillary column PERMABOND CW 20 M ID: 0.53 mm, film thickness: 0.5 µm length: 10 m	1 pz
MNL72351525	Capillary column PERMABOND CW 20 M ID: 0.53 mm, film thickness: 0.5 µm length: 25 m	1 pz
MNL72351710	Capillary column PERMABOND CW 20 M ID: 0.53 mm, film thickness: 2.0 µm length: 10 m	1 pz
MNL72351725	Capillary column PERMABOND CW 20 M ID: 0.53 mm, film thickness: 2.0 µm length: 25 m	1 pz



Code	Description	Packaging
MNL72351730	Capillary column PERMABOND CW 20 M ID: 0.53 mm, film thickness: 2.0 μ m Length: 30 m	1 pz
MNL72354910	Capillary column PERMABOND CW 20 M ID: 0.53 mm, film thickness: 1.0 μ m length: 10 m	1 pz
MNL72354925	Capillary column PERMABOND CW 20 M ID: 0.53 mm, film thickness: 1.0 μ m length: 25 m	1 pz
MNL72354930	Capillary column PERMABOND CW 20 M ID: 0.53 mm, film thickness: 1.0 μ m Length: 30 m	1 pz
MNL72382710	Capillary column PERMABOND CW 20 M ID: 0.32 mm, film thickness: 0.35 μ m length: 10 m	1 pz
MNL72382725	Capillary column PERMABOND CW 20 M ID: 0.32 mm, film thickness: 0.35 μ m length: 25 m	1 pz
MNL72382750	Capillary column PERMABOND CW 20 M ID: 0.32 mm, film thickness: 0.35 μ m length: 50 m	1 pz