



Colonne capillari per analisi solventi e alcoli PERMABOND CW 20 M

Macherey Nagel

Colonne capillari per analisi solventi e alcoli PERMABOND CW 20 M.



| Code        | Description                                                                                                | Packaging |
|-------------|------------------------------------------------------------------------------------------------------------|-----------|
| MNL72306010 | Capillary column PERMABOND CW 20 M ID: 0.25 mm, film thickness: 0.25 µm length: 10 m                       | 1 pz      |
| MNL72306025 | Capillary column PERMABOND CW 20 M ID: 0.25 mm, film thickness: 0.25 µm length: 25 m                       | 1 pz      |
| MNL72306030 | Capillary column PERMABOND CW 20 M ID: 0.25 mm, film thickness: 0.25 µm length: 30 m                       | 1 pz      |
| MNL72306050 | Capillary column PERMABOND CW 20 M ID: 0.25 mm, film thickness: 0.25 µm length: 50 m                       | 1 pz      |
| MNL72306060 | Capillary Column PERMABOND CW 20 M ID: 0.25 mm, film thickness: 0.25 µm Length: 60 m                       | 1 pz      |
| MNL72306410 | Capillary column PERMABOND CW 20 M ID: 0.10 mm, film thickness: 0.10 µm Length: 10 m                       | 1 pz      |
| MNL72329610 | Capillary column PERMABOND CW 20 M ID: 0.32 mm, film thickness: 0.5 µm length: 10 m                        | 1 pz      |
| MNL72329625 | Capillary column PERMABOND CW 20 M ID: 0.32 mm, film thickness: 0.5 µm length: 25 m                        | 1 pz      |
| MNL72329630 | Capillary column PERMABOND CW 20 M ID: 0.32 mm, film thickness: 0.5 µm length: 30 m                        | 1 pz      |
| MNL72329650 | Capillary column PERMABOND CW 20 M ID: 0.32 mm, film thickness: 0.5 µm length: 50 m                        | 1 pz      |
| MNL72329660 | Capillary column PERMABOND CW 20 ID: 0.32 mm, film thickness: 0.5 µm length: 60 m                          | 1 pz      |
| MNL72332110 | Capillary column PERMABOND CW 20 M ID: 0.32 mm, film thickness: 0.25 µm length: 10 m                       | 1 pz      |
| MNL72332125 | Capillary column PERMABOND CW 20 M ID: 0.32 mm, film thickness: 0.25 µm length: 25 m                       | 1 pz      |
| MNL72332130 | Capillary column PERMABOND CW 20 M ID: 0.32 mm, film thickness: 0.25 µm Length: 30 m                       | 1 pz      |
| MNL72332150 | Capillary column PERMABOND CW 20 M ID: 0.32 mm, film thickness: 0.25 µm length: 50 m                       | 1 pz      |
| MNL72332160 | Capillary column PERMABOND CW 20 M ID: 0.32 mm, film thickness: 0.25 µm length: 60 m Special manufacturing | 1 pz      |
| MNL72351510 | Capillary column PERMABOND CW 20 M ID: 0.53 mm, film thickness: 0.5 µm length: 10 m                        | 1 pz      |
| MNL72351525 | Capillary column PERMABOND CW 20 M ID: 0.53 mm, film thickness: 0.5 µm length: 25 m                        | 1 pz      |
| MNL72351710 | Capillary column PERMABOND CW 20 M ID: 0.53 mm, film thickness: 2.0 µm length: 10 m                        | 1 pz      |
| MNL72351725 | Capillary column PERMABOND CW 20 M ID: 0.53 mm, film thickness: 2.0 µm length: 25 m                        | 1 pz      |



| Code        | Description                                                                          | Packaging |
|-------------|--------------------------------------------------------------------------------------|-----------|
| MNL72351730 | Capillary column PERMABOND CW 20 M ID: 0.53 mm, film thickness: 2.0 µm Length: 30 m  | 1 pz      |
| MNL72354910 | Capillary column PERMABOND CW 20 M ID: 0.53 mm, film thickness: 1.0 µm length: 10 m  | 1 pz      |
| MNL72354925 | Capillary column PERMABOND CW 20 M ID: 0.53 mm, film thickness: 1.0 µm length: 25 m  | 1 pz      |
| MNL72354930 | Capillary column PERMABOND CW 20 M ID: 0.53 mm, film thickness: 1.0 µm Length: 30 m  | 1 pz      |
| MNL72382710 | Capillary column PERMABOND CW 20 M ID: 0.32 mm, film thickness: 0.35 µm length: 10 m | 1 pz      |
| MNL72382725 | Capillary column PERMABOND CW 20 M ID: 0.32 mm, film thickness: 0.35 µm length: 25 m | 1 pz      |
| MNL72382750 | Capillary column PERMABOND CW 20 M ID: 0.32 mm, film thickness: 0.35 µm length: 50 m | 1 pz      |



**CARLO ERBA Reagents S.r.l.**  
 VIA R. MERENDI, 22  
 20007 CORNAREDO (MI)  
 TEL. +39 02 93 99 190  
 FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

