



**Luminaris color Hi green Hi Rox qPCr Master Mix**

**Thermo Scientific**

Luminaris color Hi green Hi Rox qPCr Master Mix 500x20ul



| Code      | Description                                             | Packaging    |
|-----------|---------------------------------------------------------|--------------|
| FL50K0362 | Luminaris color Hi green Hi Rox qPCr Master Mix500x20ul | 500x20ul (*) |



**CARLO ERBA Reagents S.r.l.**  
 VIA R. MERENDI, 22  
 20007 CORNAREDO (MI)  
 TEL. +39 02 93 99 190  
 FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

