



Standard per ICP-MS, alluminio, 1000 µg/mL (1000 ppm) in 2-5% acido nitrico (HNO₃)

Reagecon Diagnostics

Standard per ICP-MS, alluminio, 1000 µg/mL (1000 ppm) in 2-5% acido nitrico (HNO₃).



Code	Description	Packaging
LLG04011734	Standard per ICP-MS, alluminio, 1000 µg/mL (1000 ppm) in 2-5% acido nitrico (HNO ₃)	100 ml



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

