

**Test in cuvetta CSB/COD Vario 0-15000 mg/l, 150 test UN 2922, 8, II, (E)****Tintometer**

Cuvette test CSB/COD Vario.

- Range: 0 - 15000mg/L.
- Parametro COD
- Compatibile con Lovibond® Tintometer® e colorimetri Hach®.
- Conforme ai requisiti US EPA.
- Utilizzabile in Lovibond® Tintometer® SC400, MC500 e SP600.



Code	Description	Packaging
LLG04667791	Test in cuvetta CSB/COD Vario 0-15000 mg/l, 150 test UN 2922, 8, II, (E)	150 pz



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

