

**Mascherina pieghevole uvex silv-Air premium 5210+, FFP2****Uvex**

Mascherina pieghevole uvex silv-Air premium 5210+, FFP2.

- Colore: bianco.
- Con valvola di espirazione.



| Code        | Description                      | Packaging |
|-------------|----------------------------------|-----------|
| LLG04719690 | uvex silv-Air premium 5210+ FFP2 | 15 pz     |



**CARLO ERBA Reagents S.r.l.**  
 VIA R. MERENDI, 22  
 20007 CORNAREDO (MI)  
 TEL. +39 02 93 99 190  
 FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

