

**Siringa da microlitro 1710 SL (22s/51/2)L, 100 µl****Hamilton**

Siringa a tenuta di gas, ago rimovibile lungh. 51 mm. Materiali a contatto con il liquido: SST, vetro borosilicato, PTFE.



Code	Description	Packaging
LLG06079705	Siringa da microlitro 1710 SL (22s/51/2)L, 100 µl	1 pz



CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

