



Provetta AF228/T Ø 10,65 mm per AF 228, AF 229

The Tintometer

Provetta per essiccazione AF228/T Ø 10,65 mm per AF 228, AF 229. Per cella coulometrica.



Code	Description	Packaging
LLG06085808	AF228/T Provetta Ø 10,65 mm per AF 228, AF 229	1 pz



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

