



Raccordi rapidi a pressione KECK

Duran

Raccordi rapidi a pressione KECK. In POM. Colore bianco. Parti di ricambio.



Code	Description	Packaging
LLG06202117	KECK Raccordi rapidi a pressione KK, bianco, KECK-ART.-No. 30-01	1 pz



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

