

**Portacampione per contaglobuli COULTER COUNTER®****Kartell**

Contenitori monouso in Polistirene, anche per contaglobuli Toa® e contaglobuli Royco- Hycl®.
Dotati di tappo in Polietilene.
Certificati CE-IVD.



Code	Description	Packaging
LLG06205689	Portacampione per contaglobuli COULTER COUNTER®, Ø 30 mm, h 56 mm	200 pz



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

