



Provette con tappo a vite, 30 ml, (LxØ): 84x30 mm

Sarstedt

Screw cap tube, 30 ml, (LxØ): 84x30 mm,



Code	Description	Packaging
LLG06226395	Provette con tappo a vite, 30 ml, (LxØ): 84x30 mm, fondo piatto, PP, pressione standard, con chiusura montata, chiusura sterile	450 pz



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

