

**Inserto in alluminio Dri-Block®, 12 x 16 mm, fondo conico****Stuart**

- Blocchi in alluminio.
- Dimensioni 95 x 75 x 50 mm.
- Profondità del foro: 48 mm.
- Eccellente conducibilità termica.



Code	Description	Packaging
LLG06240679	Inserto in alluminio Dri-Block®, 12 fori x prov. diam. 16 mm, fondo conico	1 pz



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

