



Termometro digitale tipo 12010/E -50...+150°C, con memoria Max/Min, incl. sonda con punta in acciaio inox 3,5x105 mm

Ludwig Schneider

Digital thermometer type 12010/E -50...+



Code	Description	Packaging
LLG06268132	Termometro digitale tipo 12010/E -50...+150°C, con memoria Max/Min, incl. sonda con punta in acciaio inox 3,5x105 mm	1 pz



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

