



Desktop mount, PE-HD, electrostatic conductive upper thread GL45(m), lower thread GL25(m), w. counter nut, silicone washer

SCAT

Desktop mount, PE-HD, electrostatic conductive upper thread GL45(m), lower thread GL25(m), w. counter nut, silicone washer.



Code	Description	Packaging
LLG06282883	Desktop mount, PE-HD, electrostatic conductive upper thread GL45(m), lower thread GL25(m), w. counter nut, silicone washer	1 pz



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

