



Test per ferro Spectroquant

Merck

Iron Test Spectroquant



Code	Description	Packaging
LLG06283155	Test del ferro Spectroquant® Metodo: fotometrico, kit per test chimici UN 3316, 9, II, (E), MR: 0,005- 5,00mg/l, 250 test	250 test



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

