



Lavavetreria da laboratorio automatica PG 8593 AE-WW-AD, acciaio inox, 835x600x600mm,
Qmax.=400 L/min

Miele

Lavavetreria da laboratorio con dosaggio liquido e asciugatura aria calda DryPlus.
Collegabile di serie tramite WiFi o Ethernet.





CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

