



Sensore di pH con sensore di temperatura SE 558X/1-NMSN, pH: 0...14, Temp. -5...100°C, lungh.: 120 mm **Knick Elektronische**

Sensore di pH con sensore di temperatura SE 558X/1-NMSN pH: 0...14, Temp. -5...100°C L: 120 mm.





CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

