



Dischi portacampioni, in alluminio, rotondo per DAB, DBS 60-3, TDIS/TDLG-A

Kern & Sohn

Piatti per campioni monouso in alluminio, Ø 90 mm.



| Code        | Description                                                                | Packaging |
|-------------|----------------------------------------------------------------------------|-----------|
| LLG06329015 | Dischi portacampioni, in alluminio, rotondo per DAB, DBS 60-3, TDIS/TDLG-A | 1 pz      |



**CARLO ERBA Reagents S.r.l.**  
 VIA R. MERENDI, 22  
 20007 CORNAREDO (MI)  
 TEL. +39 02 93 99 190  
 FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

