



Electronical Digital Thermometer ad 15 lab -50...+200:0,1°C, semiconductor sensor stainless steel,
35x5mm, 500mm PTFE cable heat resistant up

Amarell



Code	Description	Packaging
LLG06490740	Termometro digitale elettronico 15 lab	1 pz



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

