

**Siringhe Combitips advanced Box 2.0, epQuality****Eppendorf**

Siringhe Combitips advanced Box 2.0, epQuality.

Per dispensatori a spostamento positivo.

Adatta anche a liquidi viscosi e volatili.



Code	Description	Packaging
LLG06492975	Siringhe Combitips advanced Box 2.0, epQuality	1 pz



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

