



Olio di calibrazione al silicone, viscosità 50 mPa.s, volume 250 ml, temperatura 23 °C

Lamy Rheology

Olio di calibrazione al silicone, viscosità 50 mPa.s, volume 250 ml, temperatura 23 °C.



Code	Description	Packaging
LLG06522055	Olio di calibrazione al silicone, viscosità 50 mPa.s, volume 250 ml, temperatura 23 °C	1 pz



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

